

Effects of DNA Methylation on Restriction







IN LABORATORY 4, THE *EcoRI* METHYLATION SYSTEM is used to illustrate the sequence specificity of a modifying enzyme that protects DNA from restriction enzyme digestion. *EcoRI* methylase adds a methyl group to the second adenine residue in the *EcoRI* recognition site, thus preventing the endonuclease from binding and cutting the DNA. *S*-Adenosyl methionine (SAM), included in the methylation reaction, donates the methyl group that is attached to the DNA molecule by the methylase.


Three samples of bacteriophage λ DNA are incubated at 37°C with *EcoRI* methylase, one of which is subsequently incubated with *EcoRI* and another of which is incubated with *HindIII*. The third sample, a control, is incubated without a restriction enzyme. Three control samples of nonmethylated DNA are also incubated with *EcoRI*, *HindIII*, and no enzyme. All of the samples are separated by electrophoresis in an agarose gel and stained. Comparison of the band patterns reveals that the methylated DNA is protected from digestion by *EcoRI*. However, methylation at the *EcoRI* site has no effect on the activity of the restriction enzyme *HindIII*.



Equipment and materials for this laboratory are available from the Carolina Biological Supply Company (see Appendix 1).

I. Set Up Methylase Reaction


ADD


M-E-	M+E-	M-E+	M+E+	M-H+	M+H+
					
λ DNA Buf/SAM H ₂ O	λ DNA Buf/SAM H ₂ O Methyl	λ DNA Buf/SAM H ₂ O	λ DNA Buf/SAM Methyl	λ DNA Buf/SAM H ₂ O	λ DNA Buf/SAM Methyl

MIX all tubes 

INCUBATE all tubes  37°C 





II. Cast 0.8% Agarose Gel


POUR gel 



SET 

III. Set Up Restriction Reaction


ADD


M-E+	M+E+	M-H+	M+H+
			
<i>EcoRI</i>	<i>EcoRI</i>	<i>HindIII</i>	<i>HindIII</i>


MIX 


INCUBATE all tubes  37°C 

IV. Load Gel and Separate by Electrophoresis


ADD to all tubes  Loading dye


LOAD gel 


ELECTROPHORESE 100–150 volts 





V. Stain Gel and View (Photograph)

STAIN gel 

RINSE gel 

VIEW gel 

PHOTOGRAPH gel 



PRELAB NOTES

Review Prelab Notes in Laboratory 3, DNA Restriction Analysis.

S-Adenosyl Methionine

S-Adenosyl methionine (SAM) is incorporated into a 2x restriction buffer, so that the same buffer is used for *both* methylation and restriction reactions. Because SAM is not very stable, mix the buffer/SAM solution just prior to the lab and discard after use. In addition, make sure to use a fresh stock of SAM not more than several months old.

To Avoid Confusion

This laboratory has two distinct steps involving two similar-sounding reagents. In the first step, DNA is preincubated with *EcoRI* methylase. In the second step, the methylated DNA is incubated with *EcoRI* restriction enzyme. To avoid mishaps, do not set out the endonuclease *EcoRI* restriction enzyme until *after* the methylation reactions are set up.

For Further Information

The protocol presented here is based on the following published methods:

- Aaij C. and Borst P. 1972. The gel electrophoresis of DNA. *Biochim. Biophys. Acta* **269**: 192–200.
- Greene P.H., Poonian M.S., Nussbaum A.L., Tobias L., Garfin D.E., Boyer H.W., and Goodman H.M. 1975. Restriction and modification of a self-complementary octanucleotide containing the *EcoRI* substrate. *J. Mol. Biol.* **99**: 237–261.
- Helling R.B., Goodman H.M., and Boyer H.W. 1974. Analysis of R-*EcoRI* fragments of DNA from lambdoid bacteriophages and other viruses by agarose-gel electrophoresis. *J. Virol.* **14**: 1235–1244.
- Sharp P.A., Sugden B., and Sambrook J. 1973. Detection of two restriction endonuclease activities in *Haemophilus parainfluenzae* using analytical agarose–ethidium bromide electrophoresis. *Biochemistry* **12**: 3055–3063.

PRELAB PREPARATION

Before performing this Prelab Preparation, please refer to the cautions indicated on the Laboratory Materials list.

1. Prepare 2x restriction buffer plus SAM within 1–2 days before the experiment.
2. Prepare aliquots for each experiment:
 - 30 μl of 0.1 $\mu\text{g}/\mu\text{l}$ λ DNA (store on ice)
 - 40 μl of 2x restriction buffer/SAM (store on ice)
 - 5 μl of *EcoRI* methylase (store on ice)
 - 3 μl each of *EcoRI* and *HindIII* (store on ice)
 - 500 μl of distilled water
 - 500 μl of loading dye

3. Prepare 0.8% agarose solution (40–50 ml per experiment). Keep agarose liquid in a hot-water bath (~60°C) throughout lab. Cover with aluminum foil to retard evaporation.
4. Prepare 1× Tris/Borate/EDTA (TBE) buffer for electrophoresis (400–500 ml per experiment).
5. Prepare ethidium bromide or methylene-blue-staining solution (100 ml per experiment).
6. Adjust water bath to 37°C.

MATERIALS

REAGENTS

Agarose (0.8%)
 Distilled water
 Ethidium bromide▼(1 µg/ml) (or
 0.025% methylene blue▼)
 λ DNA (0.1 µg/µl)
 Loading dye
 2× Restriction buffer/SAM▼
 Enzymes:
EcoRI
HindIII
EcoRI methylase
 1× Tris▼/Borate/EDTA (TBE) buffer

SUPPLIES AND EQUIPMENT

Aluminum foil
 Beakers for agarose and for
 waste/used tips
 Camera and film (optional)
 Electrophoresis box
 Latex gloves
 Masking tape
 Microfuge (optional)
 Micropipettor (0.5–10-µl) + tips
 Parafilm or wax paper (optional)
 Permanent marker
 Plastic wrap (optional)
 Power supply
 Test tube rack
 Transilluminator (optional)▼
 Tubes (1.5-ml)
 Water baths (37°C and 60°C)

▼ See Appendix 4 for Cautions list.

METHODS

I. Set Up Methylase Reaction

(30 minutes, including incubation)



1. Use a permanent marker to label six 1.5-ml tubes, in which methylation and restriction reactions will be performed:

M– E– = no methylase, no *EcoRI*
 M+ E– = methylase, no *EcoRI*
 M– E+ = no methylase, *EcoRI*
 M+ E+ = methylase, *EcoRI*
 M– H+ = no methylase, *HindIII*
 M+ H+ = methylase, *HindIII*


M+E+	4 μ l	5 μ l	—	1 μ l
M-H+	4 μ l	5 μ l	1 μ l	—
M+H+	4 μ l	5 μ l	—	1 μ l

It is not necessary to change tips when adding the same reagent. The same tip may be used for all tubes, provided the tip has not touched the solution already in the tubes.

To avoid confusing methylase with reagents in Part III, discard after completing Step 7.



After several hours, methylase loses activity and the reaction stops.



Gel is cast directly in the box in some electrophoresis apparatuses.

Too much buffer will chan-

3. Collect and place reagents in the test tube rack on the lab bench.
4. Add 4 μ l of DNA to each reaction tube. Touch the micropipettor tip to the side of the reaction tube, as near to the bottom as possible, to create capillary action to pull the solution out of the tip.
5. Use a *fresh tip* to add 5 μ l of restriction buffer/SAM to a clean spot on each reaction tube.
6. Use a *fresh tip* to add specified volume of distilled water to appropriate tubes.
7. Use a *fresh tip* to add 1 μ l of *Eco*RI methylase to Tubes M+E⁻, M+E⁺, and M+H⁺.
8. Close tube tops. Pool and mix reagents by pulsing in a microfuge or sharply tapping the tube bottom on the lab bench.
9. Place the reaction tubes in a 37°C water bath, and incubate for a minimum of 20 minutes. Reactions can be incubated for longer periods of time.



Following incubation, freeze reactions at -20°C until ready to continue. Thaw reactions before continuing to Section III, Step 1.

II. Cast 0.8% Agarose Gel

(5 minutes)

1. Seal the ends of the gel-casting tray with tape, and insert a well-forming comb. Place gel-casting tray out of the way on the lab bench, so that agarose poured in next step can set undisturbed.
2. Carefully pour enough agarose solution into the casting tray to fill to a depth of about 5 mm. Gel should cover only about one-third the height of comb teeth. Use a pipette tip to move large bubbles or solid debris to the sides.

removed from a dry ice can
near the wells.





III. Set Up Restriction Reaction

(30 minutes, including incubation)

1. Remove methylation reactions from the water bath or thaw the tube stored in the freezer.
2. Use the matrix below as a checklist while adding reagents to each reaction. Read down each column, adding the same reagent to all appropriate tubes. Use a fresh tip for each reagent. Refer to detailed instructions that follow.

Tube	<i>Eco</i> RI	<i>Hind</i> III
M-E-	—	—
M+E-	—	—
M- E+	1 μ l	—
M+E+	1 μ l	—
M-H+	—	1 μ l
M+H+	—	1 μ l


3. Collect *Eco*RI and *Hind*III, and place them on ice on the lab bench.
4. Add 1 μ l of *Eco*RI to Tubes M-E+ and M+E+.
5. Use a fresh tip to add 1 μ l of *Hind*III to tubes labeled M-H+ and M+H+.
6. Close tube tops. Pool and mix reagents by pulsing in a microfuge or sharply tapping the tube bottom on the lab bench.
7. Place the reaction tubes in a 37°C water bath, and incubate restriction reactions for a minimum of 20 minutes. Reactions can be incubated for longer periods of time.



After several hours, enzymes
lose activity and the reaction
stops.




Following incubation, freeze reactions at -20°C until ready to continue. Thaw reactions before continuing to Section IV, Step 1.



A piece of dark construction paper beneath the gel box will give the wells more visibility.


- b. Place six individual droplets of loading dye (1 μl each) on a small square of Parafilm or wax paper. Withdraw contents from the reaction tube, and mix with a loading dye droplet by pipetting in and out. Immediately load dye mixture according to Step 2. Repeat successively, *with a clean tip*, for each reaction.
2. Use a micropipettor to load 10 μl of each reaction tube into a separate well in the gel, as shown in diagram below. *Use a fresh tip for each reaction.*
 - a. Use two hands to steady the micropipettor over the well.
 - b. Before loading sample, make sure that there are no bubbles in the wells. If bubbles exist, remove them with a micropipettor tip.
 - c. If there is air in the end of the tip, carefully depress the plunger to push the sample to the end of the tip. (If an air bubble forms a “cap” over the well, DNA/loading dye will flow into buffer around edges of well.)
 - d. Dip the micropipettor tip through the surface of the buffer, center it over the well, and gently depress the micropipettor plunger to slowly expel the sample. Sucrose in the loading dye weighs down the sample, causing it to sink to the bottom of the well. *Be careful not to punch the tip of the micropipettor through the bottom the gel.*



Alternately, set power supply on lower voltage, and run the gel for several hours. When running two gels from the same power supply, the current is double that for a single gel at the same voltage.

3. Close the top of the electrophoresis box, and connect the electrical leads to a power supply, anode to anode (red-red) and cathode to cathode (black-black). Make sure that both electrodes are connected to the same channel of power supply.
4. Turn power supply on, and set to 100–150 volts. The ammeter should register approximately 50–100 milliamperes. If current is not detected, check connections and try again.
5. Separate by electrophoresis for 40–60 minutes. Good separation will have occurred when the bromophenol blue band has moved 4–7 cm from the wells. If time allows, carry out electrophoresis until the bromophenol blue band nears the end of the gel. *Stop* electrophoresis before bromophenol blue band runs off end of gel.
6. Turn off power supply, disconnect leads from the inputs, and remove top of electrophoresis box.
7. Carefully remove casting tray from electrophoresis box, and slide the gel into a disposable weigh boat or other shallow tray. Label staining tray with your name.
8. Stain and view gel using one of the methods described in Sections VA and VB.

Staining may be performed by an instructor in a controlled area when students are not present.



Staining time increases markedly for thicker gels.

Ethidium bromide solution may be reused to stain 15 or more gels. When staining time increases markedly, dispose of ethidium bromide solution as explained in laboratory 3.

CAUTION

Review Responsible Handling of Ethidium Bromide in Laboratory 3. Wear latex gloves when staining, viewing, and photographing gels and during cleanup. Confine all staining to a restricted sink area. For further information, see Appendix 4.

1. Flood gel with ethidium bromide solution (1 $\mu\text{g/ml}$), and allow to stain for 5–10 minutes.
2. Following staining, use a funnel to decant as much ethidium bromide solution as possible from the staining tray back into the storage container.
3. Rinse the gel and tray under running tap water.
4. If desired, gels can be destained in tap water or distilled water for 5 minutes or more to remove background ethidium bromide.




Staining intensifies dramatically if rinsed gels set overnight at room temperature. Stack staining trays, and cover top gel with plastic wrap to prevent desiccation.

5. View under UV transilluminator or other UV source.

CAUTION

Ultraviolet light can damage eyes. Never look at unshielded UV light source with naked eyes. View only through a filter or safety glasses that absorb harmful wavelengths. For further information, see Appendix 4.

6. Photograph with a Polaroid or digital camera.
7. Take time for responsible cleanup.
 - a. Wipe down camera, transilluminator, and staining area.
 - b. Wash hands before leaving lab





5. View gel over light box; cover surface with plastic wrap to prevent staining.
6. Photograph with a Polaroid or digital camera.


RESULTS AND DISCUSSION

Each type II restriction enzyme has a corresponding methylase that recognizes the same nucleotide sequence. However, rather than cutting the DNA at that point, the methylase adds a methyl group within the recognition sequence. This “modification” blocks the restriction endonuclease from recognizing and binding to the restriction site. Within the bacterium, methylation protects the host DNA from cleavage by its own restriction enzyme. Unmethylated foreign DNA is not protected from cleavage.

The methylation reaction requires *S*-adenosyl methionine (SAM), the universal methyl-donating molecule in both prokaryotes and eukaryotes. As its name implies, SAM is composed of the nucleoside adenosine and the amino acid methionine. The donation of a methyl from the methionine portion of the molecule converts it into *S*-adenosyl homoserine.

One common occurrence in this laboratory is partial methylation, where methyl groups are added to only a fraction of the *EcoRI* sites within the λ DNA molecules. Cleavage at the unprotected sites produces a partial digest, yielding restriction fragments of varying lengths. These fragments are evidenced as lower-molecular-weight bands in an agarose gel. The intensity of the bands is inversely proportional to the level of DNA methylation.

1. Examine the photograph of your stained gel (or view on a light box or overhead projector). Compare your gel to the ideal gel shown on the next page. How can you account for differences in separation and band intensity?
2. What does the M+H+ control tell you about *EcoRI* methylation?
3. What does the M+E- control tell you about methylation?



Ideal Gel

Incomplete Methylation

Faint bands in Lane M+E+ DNA partially cut by *Eco*RI

the pAMP “backbone” also contains an *Eco*RI site, which is not the intended cloning site for the *Eco*RI fragment in Step 2.

- a. Draw a diagram of this cloning experiment.
 - b. Explain how *Eco*RI methylase could be used to solve this experimental problem.
6. Which nucleotide(s) is (are) methylated by *Eco*RI methylase? Draw the structure of the newly methylated base.

FOR FURTHER RESEARCH

1. Design and execute a series of experiments to study the kinetics of a methylation reaction.
 - a. Determine the approximate percentage of sites protected at various time points.
 - b. Repeat the experiments with several methylase dilutions and several DNA dilutions.
 - c. In each case, at what time point does protection appear to be complete?
2. Design and execute experiments to use *Eco*RI methylase to map the locations of *Eco*RI restriction sites in the λ genome.
3. Research the use of methylases in constructing a genomic library.
4. Research the role of DNA methylation in gene regulation in higher organisms.
5. Research the role of methylation in controlling the movement of transposable elements in maize (corn).